




AC- 609120	CUSTOMER DRAWING	Page 1 of 1
Revision 02	PS CrystalQuick™ Plate , 96 Well, Solid Bottom, Round Well, Item-No. 609 120	



Cat.-No.	Description	A	B	C	D	E	F	G	H	I	J	K
609 120	Crystal Quick™	127.76 ±0.2	85.48 ±0.2	14.4 ±0.2	Ø 2.13	2.5	11	14.38	12.84	2.99	9	1

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date January 04, 2008	Date January 04, 2008	Date January 05, 2008	
Date July 29, 2005	Name F. Falter	Name F. Sonntag	Name F. Sonntag	