




AC- 609171/ 609180	<b>CUSTOMER DRAWING</b>	Page 1 of 1
Revision 03	<b>96 Well CrystalQuick™ / CrystalQuick™Plus,PS</b> Item-No. 609 171/ 609 180	 <b>greiner bio-one</b>

NOT TO SCALE



Cat.-No.	Description	A	B	C	D	D'	E	F	G	H	I	J	K
609 171	CrystalQuick™	127.76 $\pm 0.2$	85.48 $\pm 0.2$	8.0 $\pm 0.15$	8.3	2	2.5	5.5	14.38	13.09	8.59	9	0.5
609 180	CrystalQuick™ Plus	127.76 $\pm 0.2$	85.48 $\pm 0.2$	8.0 $\pm 0.15$	8.3	2	2.5	5.5	14.38	13.09	8.59	9	0.5

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date October 13, 2008	Date October 13, 2008	Date October 14, 2008	
Date June 03, 2008	Name F. Falter	Name F. Sonntag	Name F. Sonntag	